## **MR NEURO DRAPES - LONG INSTRUCTIONS FOR USE**

C 0123

---

## Table of Contents

- I. Intended Use**
- II. Device Description**
- III. General Warnings and Precautions**
- IV. Use Instructions**
  - A. Preparation**
  - B. Hanging the Drape**
  - C. Post Procedure**

- I. **Intended Use:** The MRI Procedure Drape is intended to provide a sterile, disposable covering of the MRI instrumentation during surgical procedures conducted under MRI imaging.

**Warning:** This device is intended for single use only. Contents of unopened, undamaged package are sterile. Do not re-sterilize.

**Caution:** Federal Law restricts this device to sale by or on the order of a physician.

II. **Device Description**

**Package Contents:**

**NGS-PD-02-L: MR Neuro Drape – Long**

MR Neuro Drape – Long and Marker Pen

**NGS-PD-03-L: MR Neuro Drape w/Extension – Long**

MR Neuro Drape – Long w/Extension and Marker Pen

**NGS-PD-02-L-E: MR Neuro Drape – Long - Europe**

MR Neuro Drape – Long and Marker Pen

**NGS-PD-03-L-E: MR Neuro Drape w/Extension – Long - Europe**

MR Neuro Drape – Long w/Extension and Marker Pen

The drape is constructed with opaque cloth and clear plastic film to create a sterile barrier for a surgical procedure. The adhesive section at the end of the clear plastic film allows attachment to the patient's scalp. The configuration includes a suspension system that allows the clear film area to travel as the patient is moved in and out of the MRI scanner bore. There are two pouches and straps at the proximal end of the drape that can be used to support cables and store instruments when not in use.

III. **General Warnings and Precautions**

**Caution:** Follow standard sterile techniques when handling the drape and any included accessories.

**Caution:** Use care not to discard the other items that are included in the drape box. These may be necessary to complete the procedure.

**Caution:** For best adhesion, avoid contact with the adhesive prior to attaching to the scanner.

**Warning:** Single Use Only. Do not attempt to re-sterilize the drape under any circumstance.

**Warning:** There are no known and reliable means of cleaning, disinfecting, repairing, and sterilizing these devices that returns them to original specifications and renders them safe and effective for reuse.


**Note:** Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and appropriate Competent Authority of the Member State in which the incident occurred.

**Note:** Safe disposal of the device: The device shall be treated as biohazardous materials and shall be disposed of accordingly per hospital policy.

IV. **Use Instructions:**

A. **Preparation**

1. The manufacturer recommends four persons for proper installation of this drape.
  - Patient should be positioned and moved through the magnet such that the head is near the back end and the body is located within the bore.
  - Three persons (two sterile and one non-sterile) should be located at the head of the patient.
  - The fourth person (non-sterile) should be located at the foot of the bed on the opposite side of the bore.
2. Five hooks must be pre-mounted on each end of the Magnetic Resonance Imaging (MRI) Magnet as in diagram below. It is recommended that the hooks be placed 12-18 inches from the nearest edge of the scanner bore. **See Figure 1.**

**Figure 1: Hook Locations**

3. Facilities using a GE Scanner, with a "Dog House" over the drive mechanism for the scanner bed, shall place a ClearPoint Neuro GE Housing Cover over this mechanism's housing prior to attaching the Drape. Ensure that the side marked "Coil Connector Side" faces the opening to the scanner. **See Figure 2.**

**Figure 2**

4. Ensure that the cover is fully seated. It should not rock under mild pressure. **See Figure 3.**

**Figure 3**

Remove the Marker Pen in the zip lock bag that is attached to the Drape pouch and set aside. The Marker Pen is in a sterile barrier pouch. The contents of the pouch are sterile. The outside surface of the pouch is non-sterile. The contents include one (1) Sterile Marking Pen.

#### **B. Hanging the Drape**

##### **MR Drape Instructions: Siemens and Philips Scanners**

1. Open the single sterile barrier drape pouch and remove the drape using standard hospital sterile techniques.
2. Unfold the protective outer drape to access the procedure drape.
3. Have two sterile people hold each side of the folded drape with the large adhesive portion (incise tape) facing the patient's head. **See Figure 4.**

**Figure 4: Adhesive area of drape with orientation arrow**

4. The large adhesive area is held vertically, facing the patient. Orient the drape with the #3 center tab (with elastic strap) on top or with the arrow pointing toward the ceiling (if applicable).
5. If it is necessary to cut a hole or holes through the adhesive area to provide room to access the patient's scalp at the planned entry point(s), this may be performed after attaching the adhesive area to the patient.
  - 5.1 Cut the adhesive area as needed for access to the scalp using an appropriate MR compatible device (e.g. Ti scissors and mosquito forceps). Discard any excess material.
6. Remove the adhesive backing and affix to the previously applied surgical incise tape on the patient, orienting any holes appropriately. **NOTE:** If the hole orientation is not correct, remove additional material. If any bone (mounting) screws will pass through the incise tape, a small incision should be made through the incise tape prior to screwing the base on to prevent twisting the incise tape.
  - 6.1 Allow the drape to unfold sufficiently to expose the 2 sets of five (5) tabs with attached elastic straps.
  - 6.2 If the drape has one set of black tabs and one set of white tabs, the black tabs are for the foot end of the table and the white tabs are for the patient (head) end of the table. **See Figure 5.**

**Figure 5: Black/White set of tabs for drape mounting**

- 6.3 If the drape has two sets of white tabs, the square tabs are for the foot end of the table and the circle tabs are for the patient (head) end of the table. **See Figure 6.**

**Figure 6: White set of tabs for drape mounting**

7. A sterile person will need to access the straps.
8. Tighten the knots of all elastic straps (where it attaches to the tab).
9. Have a non-sterile person positioned at each end of the scanner bore.
10. Starting with the bottom right (#1) tab, have the sterile person pass the tab to the non-sterile person at the head of the bed. This non-sterile person will then pass the tab through the bore to the non-sterile person at the foot of the bed. The person at the foot of the bed then attaches the tab in a clockwise direction to the appropriate hook that is attached to the MR magnet face. Repeat this procedure for all five (5) foot end tabs, making sure each tab is passed separately through the bore.
11. The second set of five (5) patient (head) end tabs, with attached elastic straps, are now ready to be exposed.
12. Repeating the steps used for the foot end of the scanner, the sterile person hands the tabs to the non-sterile person who attaches the tab to the appropriate hook that is attached to the MR magnet face. When completed, all elastic straps should travel straight through the bore without crossing other straps.
13. Two sterile people, one at each corner of the drape, lift up the top corners of the drape exposing the sterile side. They remove the backing of the adhesive strips along the edges of the drape, starting in the middle and working out.
14. A Sterile person affixes the drape to the patient (head) side of the scanner. Center of the drape should be affixed to the scanner above the center hook. Each hook is covered by the drape when mounted.
15. Separate the four pieces of tape with arrows on each (if present) from the clear side of the drape so it can extend its full length into the scanner.  
**Caution:** Do not tear the clear portion of the drape material while separating the pieces of tape.
16. Move patient to center of scanner bore to extend drape and check for proper functioning of elastic straps.

**MR Drape w/ Extension Instructions: Additional steps required for GE scanners.**

17. Verify the GE Housing Cover has been placed over the GE drive mechanism.
18. Unfold the protective outer drape to access the procedure drape. Orient the drape with the adhesive backing with arrow pointing up facing the patient and the folded portion that covers the table on the bottom.
19. Remove the adhesive backing and attach to the patient as previously described.
20. Un-fold the table covering portion of the drape, covering the length of the table extension and beyond the Table Drive Cover (shown below with the Table moved into the scanner bore). Do not pull the sides down until after attaching the drape to the face of the bore.
21. Attach the tabs to both ends of the scanner and attach the drape to the face of the scanner as previously described.
22. Un-fold the remaining sides of the drape on the table so they cover the sides of the table.  
**See Figure 7.**

**Figure 7: Assembled GE Drape**

**Completed Drape Assembly**

23. Perform intended procedure.

**C. Post Procedure**

1. When procedure is complete, pull adhesive edges of the drape off the magnet and roll drape up slightly to prevent any blood/fluid from leaking onto the magnet bore.
2. Remove elastic straps first from head-side hooks.
3. The sterile person cuts the drape off the patient's head.
4. Move patient out of the magnet.
5. Remove the elastic straps from foot-side hooks.
6. Carefully fold and dispose of the Drape in an appropriate method.
7. The GE Housing Cover can be cleaned with normal household or hospital disinfectants.

**Caution: The Housing Cover shall not be cleaned with alcohol.**

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ClearPoint Neuro, Inc.