# CLEARPOINT

## NEURO

### **SmartTwist MR HAND DRILL**

### **SmartTip MR DRILL BIT KITS**

This product uses Surgibit® technology

### **INSTRUCTIONS FOR USE**

*Surgibit* is a registered trademark of Surgibit IP Holdings Pty Limited and is used with permission.

## Table of Contents

|      |                                                                      |
|------|----------------------------------------------------------------------|
| I.   | <b>Intended Purpose</b>                                              |
| II.  | <b>Device Description</b>                                            |
| III. | <b>General Warnings and Precautions</b>                              |
| IV.  | <b>MRI Safety Information</b>                                        |
| V.   | <b>Use Instructions</b>                                              |
| A.   | <b>Preparation</b>                                                   |
| B.   | <b>Operating the Hand Drill</b>                                      |
| C.   | <b>Inserting the Drill Bit into the Drill</b>                        |
| D.   | <b>Removing the Drill Bit from the Hand Drill</b>                    |
| E.   | <b>Using the Drill Bit and Accessories</b>                           |
| F.   | <b>Storage and Technical Specifications</b>                          |
| VI.  | <b>Summary of Safety and Clinical Performance (SSCP) Information</b> |

### I. Intended Purpose

The SmartTwist MR Hand Drill and accessories are intended to provide access through the skull for ventriculostomy or other neurological procedures such as biopsy or catheter placement, in or near an MR scanner of 3T maximum field strength. The SmartTwist MR Hand Drill and accessories are intended to be used only when the scanner is not performing a scan. The SmartTwist MR Hand Drill is intended for single use only.

#### Contraindications

There are no contraindications for the SmartTwist MR Hand Drill and SmartTip Drill Bit.

#### Intended Users

Physicians are the only intended users of this device. Any physician trained in surgery with the use of manual cranial drills may use this device.

**Warning:** This device is intended for single use only and is Ethylene Oxide sterilized. Contents of unopened, undamaged package are sterile. Do not re-sterilize. Do not use if package is damaged or opened.

**Caution:** Federal (U.S.) law restricts this device to sale by or on the order of a physician.

### II. Device Description

#### Package Contents:

**REF** NGS-HD-01: SmartTwist MR Cranial Hand Drill  
MR Compatible Hand Drill

#### Associated Devices:

**REF** NGS-DB-45: SmartTip MR Drill Bit Kit, 4.5-mm  
4.5-mm Surgibit® Drill Bit, 3.4-mm Surgibit® Drill Bit, Drill Stop (2), Ruler, Lancet

**REF** NGS-DB-45-S: SmartTip MRI Short Drill Bit Kit, 4.5-mm  
4.5-mm Surgibit® Drill Bit, 3.4-mm Surgibit® Drill Bit, Drill Stop (2), Ruler, Lancet

**REF** NGS-DB-45-E: SmartTip MR Drill Bit Kit, 4.5-mm  
4.5-mm Surgibit® Drill Bit, 3.4-mm Surgibit® Drill Bit, Drill Stop (2), Ruler

**NOTE:** NGS-DB-45 and NGS-DB-45-S are not CE marked.

### III. General Warnings and Precautions

|                          |                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Warning statement</b> | Indicates information regarding possible injury, death, or other serious consequences which could occur with use of the device. |
| <b>Caution statement</b> | Indicates information regarding possible minor harm or damage which could occur with use of the device.                         |

**Warning:** The device is intended for single-use-only and is provided sterile. Do not re-sterilize.

**Warning:** Do not use if the package or any of the parts of the device are damaged.

**Warning:** DO NOT perform an MRI Scan while the Drill Bit or Lancet are being held or are in contact with tissue. Use the Drill

Bits or Lancet only when the MRI scanner is not running a scan. The Drill Bit and Lancet may cause injuries from RF-induced heating if they are near tissue during a scan.

**Warning:** There are no known and reliable means of cleaning, disinfecting, repairing, and sterilizing these devices that returns them to original specifications and renders them safe and effective for reuse. Infection or performance degradation may occur if device is reused.

**Warning:** Do not use the Drill Bits at speeds greater than 1500 RPM.

**Caution:** Grasp the Drill Bits only on the smooth shank section. Do not grasp the Drill Bits by the flutes or the point.

**Caution:** Do not use Drill Bits less than 2.0-mm in diameter. The Hand Drill may not retain the Bits as desired during the drilling procedure.

**Caution:** There is an inherent risk in creating access through the skull that the device may cause a subdural bleed. Monitor drilling progress closely during usage.

**Precaution:** Handle all components using standard hospital sterile practices.

**Note:** Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### IV. MRI Safety Information

- Non-clinical testing has demonstrated the MR Drill Kits are MR Conditional. They provide the clinical benefit of use in or near an MRI scanner safely under the following conditions:
  - Static Magnetic field of 1.5 Tesla and 3 Tesla
  - Spatial gradient field of 517 Gauss/cm (5.17 T/m) or less
- DO NOT perform an MRI Scan while the Drill Bits, Depth Stop, or Lancet are in the scanner. DO NOT perform an MRI Scan while the Drill Bit or Lancet are being held or are in contact with tissue. Use the Drill Bits or Lancet only when the MRI scanner is not running a scan. The Drill Bit and lancet may cause injuries from RF-induced heating if they are near tissue during a scan. REMOVE the Drill Bit, Lancet, and Depth Stop from the scanner before performing an MRI scan.

#### V. Use Instructions

##### A. Preparation

Do not use the Hand Drill if sterilization is compromised. Examine the Hand Drill for any apparent damage. Do not use the Hand Drill if it is damaged.

Turn the Drill at the expected procedure speed and confirm that the Drill Chuck is spinning freely.

Do not use the Drill Bits if sterilization is compromised. Examine the Drill Bits for any signs of damage or corrosion. Do not use a Drill Bit if it is damaged or if corrosion is present.

Do not use the Lancet if sterilization is compromised. Examine the Lancet for any signs of damage or corrosion. Do not use the Lancet if it is damaged or if corrosion is present.

##### B. Operating the Hand Drill

The MR Hand Drill operates by placing a Drill Bit in the Drill Chuck, tightening the Chuck, and then turning the Handle. These parts of the Drill are shown in **Figure 1**.


Figure 1


Figure 2


Figure 3

The Hand Drill Chuck spins as the Handle is turned. Therefore, the Handle speed directly controls the Chuck speed. To operate the Hand Drill at a higher speed, turn the Handle faster. To operate the Hand Drill at a lower speed, turn the Handle slower. The drill operates at a 3:1 gear ratio.

To turn the Drill in the forward direction, grasp the Drill and turn the Handle clockwise, as shown in **Figure 2**.

The Hand Drill may also be run in the reverse direction by turning the Handle in the counterclockwise direction, as shown in **Figure 3**.


**Caution:** Keep hands away from the inner part of the Handle. Hand or glove may be pinched.

**Caution:** Keep hands away from the bottom of the Handle Knob. Hand or glove may be pinched.


#### C. Inserting the Drill Bit into the Drill

Hold on to the Hand Drill so that the Handle cannot turn.


Hold the Handle in place and turn the Drill Chuck to the left to open it, as shown in **Figure 4**.

**Figure 4**

**Figure 5**

**Figure 6**

Turn the Drill Chuck until the collet is open far enough to accept the Drill Bit.

Insert the Drill Bit until it is fully seated inside the Chuck. Press down on the round shank section of the Drill Bit firmly to make sure it is fully seated.

**Warning:** A slight pull or twist caused by the scanner's magnetic field may be felt when handling the Drill Bit, Lancet, and Depth Stop.

**Caution:** The drill bit should not be removed from or assembled to the drill on or near the patient table due to risk of injury.

**Caution:** Grasp the Drill Bits only on the smooth shank section. Do not grasp the Drill Bits by the flutes or the point, as a cut of a glove or hand may occur.

To tighten the Chuck around the Drill Bit, hold the handle in place and turn the Chuck to the right, as shown in **Figure 5**.

If using the Drill Bit with a Drill other than the MR Hand Drill:

Make sure the Drill can accept the Drill Bit prior to the procedure.

Open the Drill's chuck per the manufacturer's instructions.

Insert the MR Drill Bit into the Chuck.

Tighten the Drill chuck per the manufacturer's instructions.

Perform a security check on the Drill Bit by grasping the round shank and pulling on it to make sure it is held securely in the Drill chuck. Drill Bit may slip during drilling if it is not secure.

#### D. Removing the Drill Bit from the Hand Drill

While holding the Handle in place, turn the Drill Chuck to the left, as shown in **Figure 6**.

When the Drill Chuck is open far enough, pull on the Drill Bit on the round shank section to remove it from the Drill.

**Caution:** Grasp the Drill Bits only on the smooth shank section. Do not grasp the Drill Bits by the flutes or the point, as a cut to the glove or hand may occur.

If using the Drill Bit with a Drill other than the MR Hand Drill:

Follow the Drill manufacturer's instructions for opening the Drill chuck.

Remove the MR Drill Bit by pulling on the round shank.

**Caution:** Grasp the Drill Bits only on the smooth shank section. Do not grasp the Drill Bits by the flutes or the point, as a cut to the glove or hand may occur.

#### E. Using the Drill Bit and Accessories

When the desired position of the hole is determined, and the Drill Bit is inserted into the Drill, place the Drill Bit point on the desired position.

Turn the Handle in the forward direction to create the hole. A light downward force may be applied during drilling.

**Warning:** A slight pull or twist caused by the scanner's magnetic field may be felt when handling the Drill Bit, Lancet, and Depth Stop.

**Warning:** DO NOT perform an MR Scan when the Drill Bit is in contact with tissue and/or being held by hand. The Drill Bit may cause the tissue to heat, resulting in injury.

Turning the Handle faster will increase the speed of the Drill Bit.


When the desired depth is reached, reverse the Handle direction and apply a light pull force to the Drill to extract the Drill Bit from the hole.

If using the MR Drill Bit with a Drill other than the MR Hand Drill, follow the Drill manufacturer's instructions for operating the Drill.

**Warning:** Do not use the MR Drill Bits at speeds greater than 1500 RPM.

If the drilling depth is determined ahead of time, the Ruler may be used to measure the distance from the tip to the desired depth on the Drill.

The Depth Stop may be attached to the Drill Bit so that the distal surface of the Drill Stop corresponds to the pre-determined depth measured in the previous step. This is shown in **Figure 7**.

**Figure 7**

Move the depth stop to the desired depth and tighten the screw until it is secure on the Drill Bit.

**Note:** If using the 4.5-mm Drill Kit, there are two depth stops supplied. The larger one is for the Drill Bit, and the smaller one is for the Lancet.

Perform a security check on the Depth Stop by pushing on it and confirming that it does not move.

**Warning:** The drill bit may go deeper than intended if the Depth Stop slips, resulting in possible patient injury.

The Depth Stop may be removed by loosening the screw and sliding the Depth Stop over the Drill Bit.

**Note:** A Lancet is provided only with the NGS-DB-45 and NGS-DB-45-S Drill Bit Kits. The Lancet may be used to pierce the dura and/or pia after the access hole is created.

**Warning:** DO NOT perform an MRI Scan when the Lancet is in contact with tissue or being held by hand. The Lancet may cause the tissue to heat, resulting in injury.

**Caution:** Dispose of the Drill, Drill Bit, Depth Stop, Lancet and Ruler in a manner consistent with best practices of the institution where the device is used. Particularly, dispose of Drill Bits in a manner consistent with best practices of sharps disposal.

#### F. Storage and Technical Specifications

##### Storage

- Store in a cool dry place

##### Technical Specifications

| Product Name                              | Length<br>(in/mm) | Diameter<br>(in/mm) | Mass<br>(g) |
|-------------------------------------------|-------------------|---------------------|-------------|
| MR Hand Drill                             | 7.5 / 191         | 3.0 / 76            | 165         |
| MR Drill Bit <sup>1</sup> , 3.4-mm        | 8.66 / 220        | 0.134 / 3.4         | 15          |
| MR Drill Bit <sup>1</sup> , 4.5-mm        | 8.66 / 220        | 0.177 / 4.5         | 25          |
| MRI Short Drill Bit <sup>1</sup> , 3.4-mm | 6.75/171          | 0.134 / 3.4         | 11          |

| Product Name                              | Length (in/mm) | Diameter (in/mm) | Mass (g) |
|-------------------------------------------|----------------|------------------|----------|
| MRI Short Drill Bit <sup>1</sup> , 4.5-mm | 6.75/171       | 0.177 / 4.5      | 18       |
| Lancet <sup>2</sup>                       | 8.66 / 220     | 0.08 / 2.0       | 12       |

<sup>1</sup>This product uses Surgibit® technology. The Surgibit® technology is protected by the following patents:

US Design Patent No. D523313, US Design Patent No. D523398, US Utilities Patent No. US 7,892,235, US Utilities Patent No. US 9,162,945, US Utilities Patent No. US 8,172,845

Lancet is not included in the NGS-DB-45-E Drill Kit

## VI. Summary of Safety and Clinical Performance (SSCP) Information

The SSCP will be made available on EUDAMED. <https://ec.europa.eu/tools/eudamed>

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ClearPoint Neuro, Inc.