# CLEARPOINT

## NEURO

### **Inflexion Head Fixation Frame - BR**

#### **INSTRUCTIONS FOR USE**

See the accompanying scanner table specific IFU for the frame provided which includes the contents of the assembly, associated devices, and instructions for attachment to the scanner table.

## Table of Contents

1. Indications for Use
2. Device Description
3. General Warnings and Precautions
4. Use Instructions
  - A. Preparation
  - B. Installation of the Fixation Screws and Base Adjustment
  - C. Positioning the Patient
  - D. Removal of Patient
  - E. Removal from MRI Unit
  - F. Head Fixation Frame Cleaning

### 1. Indications for Use

The Inflexion Head Fixation Frame is intended for use as a device to clamp and hold the patient's head in a particular position for procedures requiring Magnetic Resonance Imaging (MRI) of the brain structure and targets.

### 2. Device Description

**Warning:** This device is intended for multiple uses. Contents are non-sterile.

**Caution:** Federal (U.S.) law restricts this device to sale by or on the order of a physician.

**Note:** The Head Fixation Frame is **MR Conditional**

The Inflexion Head Fixation Frame is designed to immobilize the patient's head during the surgical procedure. The patient's head is rigidly attached to the Head Fixation Frame by the physician in the desired position (lateral, oblique, prone or supine) for the planned surgical procedure. Positioning of the head may be performed remotely with the Arc separated from the Base or with the Arc attached to the Base.

See the accompanying scanner table specific IFU for the frame provided which includes the contents of the assembly, associated devices, and instructions for attachment to the scanner table.

### 3. General Warnings and Precautions

The Head Fixation Frame (HFF) is intended for multi-use and is provided non-sterile. It should be cleaned and disinfected (e.g. using disinfectant wipes) before initial use and after each use and used in conjunction with an ClearPoint Neuro Sterile Drape. The Head Fixation Frame has been tested for 100 uses.

The following are warnings, cautions, and safety requirements that apply to the Inflexion HFF; consult the device specific instructions for all devices used in conjunction with the Inflexion HFF for warnings specific to those devices.

**Warning:** Do not clamp/affix the Head Fixation Frame to damaged or diseased bone. Only attach to stable bone to insure proper fixation.

**Warning:** Do not clamp/affix the Head Fixation Frame to patients with skull prostheses or previous reconstructive surgeries.

**Warning:** Do not over tighten the head fixation screws; select the proper clamping pressure related to the skull thickness and bone structure. Clamping to the skull only to be performed by surgeons that are trained in this aspect of neurosurgery.

**Warning:** Make sure that the patient's eyes do not contact the skull pins during head positioning. Eye contact with the pins can lead to serious injury.

**Warning:** Prior brain targeting images may not be accurate if the patient is moved from ISO center or if pressure from 10 lb<sub>r</sub> up to 25 lb<sub>r</sub> is applied to the patient's skull (which may occur when installing bone screws into the skull). Re-scan the patient to confirm the target's location with new images prior to proceeding. The following reminder can be found etched onto the face of the Lower Ring: **WARNING – RE-SCAN AFTER APPLYING FORCE**

- Caution:** This device is designed to use DORO Skull Pins. Failure to use the recommend Skull Pins may result in improper function of the Head Fixation Frame. (Recommended: DORO® Skull Pins, Manufactured by PMI - Pro Med Instruments)
- Caution:** This device is to be used only by physicians trained by ClearPoint Neuro personnel.
- Caution:** Verify all knobs and screws are secure prior to performing imaging since movement of the skull could occur during the procedure.
- Caution:** Ensure that the patient is supported with pads/foam/rolled towels etc. to support their body from putting undue pressure on the neck and skull. The neck should be in as much of a neutral position as can be achieved based on anatomy and patient positioning.
- Caution:** The Aluminum Screw Tool, the Aluminum Screw Tool Crossbar and the Titanium Hex Knob components of the Head Fixation Kit must be removed from the scanner bore before MRI scanning.
- Caution:** The compatibility of other Coil devices should be evaluated before use with the Inflexion Head Fixation Frame.
- Caution:** The Head Fixation Frame is designed to be robust; however, any attached Flex Coil (optional) is sensitive electronic equipment that requires handling with care.
- Caution:** All other Coils must be removed from the MRI table, prior to installing the ClearPoint Neuro Head Fixation Frame.
- Caution:** Notify ClearPoint Neuro immediately if any screws break or holes are stripped out on the Head Fixation Frame.
- Note:** Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and appropriate Competent Authority of the Member State in which the incident occurred.
- Note:** Safe disposal of the device: The Titanium Skull Pins shall be treated as biohazardous materials and shall be disposed of accordingly. The Head Fixation Frame shall be disposed of per hospital policy.

#### General Precautions

Clean the Head Fixation Frame (e.g. using disinfectant wipes) before initial use and after each use.

#### 1.5T & 3T Environment Compatibility (MR conditionality)

Non-clinical testing has demonstrated the Inflexion HFF is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 or 3 Tesla
- Maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of "Normal Operating Mode" scanning.

In non-clinical testing, the Inflexion Head Fixation Frame produced a temperature rise of less than 1°C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of MR scanning in a 1.5T Philips Achieva MR scanner and in a 3T Siemens Verio MR scanner.

## 4. Use Instructions

### A Preparation

The Head Fixation Frame (HFF) is intended for multi-use and is provided non-sterile. It should be cleaned and disinfected (e.g., using disinfectant wipes) before initial use and after each use and used in conjunction with an ClearPoint Neuro Sterile Drape. The Head Fixation Frame uses Head Fixation Screws and Skull Pins for patient fixation. Head fixation is accomplished using an arrangement of four to six post/fixation screw/skull pin combinations. The pins are tightened/loosened by turning the fixation screws mounted in threaded holes in the posts, using the screw tool. The Skull Pins required to fixate the patient are single-use, sterile components, manufactured by: Pro Med Instruments

- DORO® Skull Pins Titanium Tip: Item # 3006-20
- DORO® Skull Pins Titanium Tip Children: Item # 3006-30

**Warning:** Do not use the Head Fixation Frame if any of its components are damaged.

1. Inspect the Head Fixation Frame along with the (optional) Flex Coil for defects. Do not proceed if defected items cannot be replaced.
2. The following terms will be used throughout this IFU to describe the various groups of Inflexion Head Fixation Frame parts;
  - Table Mounts: The side pieces that clamp onto the scanner table and secure the assembly to the table.
  - Base (upper and lower): That portion that secures the Ring, allowing for tipping and rotation of the Ring.
  - Ring (upper and lower): The portion that holds the Posts and Fixation Screws.
3. Remove all accessory MRI Coils from the MRI table (i.e. Spine Coil, Head Coil, and others).

### B. Installation of the Fixation Screws and Base Adjustment

Recommend Skull Pins:

#### NGS-SP-01: Skull Pin Tips

MR Compatible Skull Pin Tips

(DORO Skull Pins Titanium Tip (DORO Item # 3006-20), Manufacturer: Pro Med Instruments)

#### NGS-SP-02: Pediatric Skull Pin Tips

MR Compatible Skull Pin Tips, Pediatric

(DORO Skull Pins Titanium Tip Children (DORO Item # 3006-30), Manufacturer: Pro Med Instruments)

**Warning:** Make sure that the patient's eyes do not contact the skull pins during head positioning. Eye contact with the pins can lead to serious injury.

**Note:** For additional patient safety, eye protection during positioning may be utilized. Safety eye wear is available in both adult and pediatric sizes.

**Note:** Handle the Skull Pins by the protection tube while removing from the package and placing into the Fixation Screws.

**Note:** Ensure each Skull Pin is fully seated.

**Note:** The installation of the Fixation Screws is common to all ClearPoint Neuro Inflexion Head Fixation Frames. The steps for mounting the HFF onto the patient's skull can be performed with the Ring attached to the Base assembly or separated. The steps described below are performed with the Ring attached to the Base.

1. Adjust the Base to aid fitting the patient into the frame.
  - a. Determine if the (upper) swivel part of the base needs to be raised or lowered. There are three positions. Determine the best position and insert the Green or Black Knobs into either side as shown below. Tighten until snug. (**Figure 1**)



**Figure 1**



- b. Insert the Ring into the Base and secure with the two (2) Yellow Knobs. As the patient is brought into the Ring (mounted) or when the patient (with Ring mounted to their skull) is brought to the Base, loosen the Green or Black knobs and allow the upper base to swivel to accept the patient / Ring. The Base can be swivel in a  $\pm 19^\circ$  arc.
2. The Inflexion HFF comes with 6 posts. Use 4 to 6 posts for head fixation (4 posts is the minimum number of posts that should be used). Place the posts so that the upper and lower screws provide counter balancing forces to one another. The more opposing the upper and lower screws are to one another, the more holding force will be applied.
  - The 2 acceptable post placements pictures shown below are the only acceptable configurations when not using the upper ring.
  - The unacceptable post placement pictures below demonstrate posts that are not providing opposing forces to one another.

See **Figure 2a** and **Figure 2b** below for examples.

**ACCEPTABLE**


**UNACCEPTABLE**


**Figures 2a**

3. An optional Upper Ring is available when the use of 4 Posts on the Lower Ring is not sufficient. It interlocks with the lower ring and is secured with two (2) Post Screws. The image below shows the Optional Upper Ring mounted to the Lower Ring.
  - When using the upper ring a minimum of 5 posts should be used with at least 1 empty space (slot) between each post.

**Figure 2b**

4. If using a Flex Coil (Siemens Coil shown for illustration), place the Flex Coil over the Head Fixation Screws using the openings in the Flex Coil. (**Figure 3**)

**Figure 3**

### **C. Positioning the Patient**

1. Carefully support the patient's skull in the Ring. Insert the Screw Tool through the desired hole in the Post (hex end first) until the Screw Tool touches the patient's head (**Figure 4**).


Figure 4

2. Determine the minimum length of each Fixation Screw Based on Screw Tool graduation seen at the outer face of the Post. Fixation Screws that will be near the Frame Base and those that might interfere with a flex Coil (when it is mounted on the outside of the Posts) should be a short as possible to reduce interference. (Figure 5)


Figure 5

3. While maintaining sterility and protecting the sleeves on each Skull Pin, press Skull Pins into the appropriately sized Fixation Screws.
4. Carefully thread each Fixation Screw/Skull Pin assembly through the appropriate hole in the Posts, using fingers. Follow with the Screw Tool until the pins protrude slightly from the patient side of the Posts.  
**Caution:** The tip of the Skull Pin is very sharp. Take care to protect the patient and clinician from injury once the plastic cover is removed from the skull pin.
5. Orient the flex Coil if it is going to be placed inside the Ring so that the openings on the Coil provide clearance for the Fixation Screws.
6. Carefully remove the sterile-protectant sleeves from the pins.  
**Warning:** Do not use the Screw Tool - Crossbar when tightening the pins. This may cause excessive force and may result in damage to the Head Fixation components or patient injury.
7. With an assistant holding the patient's head, gently and carefully engage the superior side of the patient's skull with the two lower skull pins. (Figure 6)


Figure 6

8. Measure the upper screw length again if necessary and engage the upper fixation screws with the patient's skull. (Figure 7)


Figure 7

9. Tighten each Fixation Screw using the Screw Tool by inserting the hex end of the Screw Tool into the Fixation Screw and rotating clockwise. The surgeon is responsible for determining the pressure required for adequate fixation.  
**Warning:** Do not over tighten the head fixation screws. Fixation of the skull should only be performed by surgeons that are trained in this aspect of neurosurgery.  
**Caution:** Ensure that the patient is supported with pads/foam/rolled towels etc. to support their body from putting undue pressure on the neck and skull. The neck should be in as much of a neutral position as can be achieved based on anatomy and patient positioning.  
**Caution:** When tightening Fixation Screws, stop turning when the "STOP" line on the Screw Tool reaches the outer face of the Post. (Figure 8) This is the minimum recommended engagement of the Fixation Screw in the Post. If more screw travel is required, remove the current Fixation Screw and replace with a longer Fixation Screw and a new Skull Pin."

Figure 8

- Note:** If the clearance between the end of the Fixation Screw and the table or other components is too small to use the Screw Tool, use the Hex Knob to tighten the screws.
10. **Note:** When a flex coil(s) is oriented on the outside of the Head Fixation Frame, longer Fixation Screws may interfere with positioning the Coil close to the frame. The Fixation Screws can be exchanged with shorter screws if necessary. Always replace the Skull Pin with a new sterile pin.
  11. Carefully inspect patient fixation in relation to the intended treatment entry point(s).
  12. If the Ring is not mounted into the Base, place it into the Base at this time.
    - a. Determine if the patient's head needs to be rotated for correct positioning for the neurosurgical procedure. If so, carefully loosen the two (2) Green or Black Knobs that secure the Upper Base.
    - b. With an assistant coordinating movement of the patient's torso, slowly rotate the patient's head and body (with the Ring rotating in the Cradle) to the desired position. Support the patient's body with padding as necessary and engage the two (2) Green or Black Knobs to secure the base and the two (2) Yellow Knobs to lock the Ring in place.

**Caution:** The patient's head must remain immobile throughout the procedure.

#### D. Removal of Patient

**Note:** Use the Screw Tool – Crossbar in the Screw Tool if additional force is required to loosen the Head Fixation Screws. Use the Hex Knob in locations that have minimum clearance.

**Caution:** Do not loosen the Upper Ring (if in use) until all Fixation Screws have been loosened. The Ring may be separated from the Base prior to loosening the screws.

1. Loosen the Fixation Screws. Normally this will be done in the reverse order that they were installed.
2. Remove the patient from Head Fixation Frame.

#### E. Removal from MRI Unit

1. Remove the Skull Pins from the Head Fixation Frame and properly dispose.
2. Loosen the top two (2) Thumbscrews that hold the Adjustable Toe Inserts.
3. Turn the front two (2) Thumbscrews of the Adjustable Toe Inserts counter clockwise until the Adjustable Toe Inserts are all the way back and then lift the Frame out of the Cradle from the adjustable toe insert side.

#### F. Head Fixation Frame Cleaning

The Head Fixation Frame (HFF) is intended for multi-use and is provided non-sterile. It should be cleaned and disinfected (e.g., using disinfectant wipes) before initial use and after each.

**Caution:** Do not use any disinfection or sterilization methods involving heat. This will damage the device. The components of the Head Fixation Frame have only been tested with cold liquid disinfectants.

1. Remove the Lower Ring, Upper Ring (*optional*), Upper Base (portion that holds the Ring), Posts, Fixation Screws and all screws for disinfection.
2. Dry the Head Fixation Frame with a lint free towel or equivalent.
3. Reassemble or place with the other components for future use. Examine the parts when cleaning and assembling to look for wear. Also check after assembly that there is no unexpected movement in the system.

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ClearPoint Neuro, Inc