Depth Stop .085 I.D.

INSTRUCTIONS FOR USE

Table of Contents

- I. Intended Use**
- II. Device Description**
- III. General Warnings and Precautions**
- IV. Use Instructions**

I. Intended Use

The MR Compatible Device Depth Stop is intended as an accessory for medical devices that may need depth limits. The device is not intended for implant. This device is intended for “single patient use only”.

Warning: This device is intended for “single patient use only”. Contents of unopened, undamaged package are sterile. Do not re-sterilize.

II. Device Description**REF NGS-DS-06: Depth Stop (6pk)**

The Depth Stop consists of two components: a ring and a set screw, both made from plastic. The package contains 6 total Depth Stops.

III. General Warnings and Precautions

Warning: Do not use if the package or any of the parts of the device are damaged.

Warning: There are no known and reliable means of cleaning, disinfecting, and sterilizing these devices that returns them to original specifications and renders them safe and effective for reuse.

Caution: The compatibility of instruments/devices with the Device Depth Stop should be evaluated before use. The components are MR safe.

Caution: Handle all components using standard hospital sterile practices. Store in cool dry place.

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IV. Use Instructions

1. Remove the device from the package. The Depth Stop is packaged in a double sterile barrier; two sealed mylar/Tyvek pouches.

2. Place the Depth Stop onto the device/instrument.

3. Carefully tighten the Depth Stop to prevent damage to the other device/instrument.

Note: Safe disposal of the device: The device shall be treated as biohazardous materials and shall be disposed of accordingly per hospital policy.

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ClearPoint Neuro, Inc.