CLEARPOINT®

NEURO

Multi-Positional Head Fixation Frame

INSTRUCTIONS FOR USE [EN]

See the accompanying scanner table specific IFU for the frame provided which includes the contents of the assembly, associated devices, and instructions for attachment to the scanner table.

Table of Contents

1. **Intended Use**
2. **Device Description**
3. **General Warnings and Precautions**
4. **Use Instructions**
 - A. **Preparation**
 - B. **Installation of the Lower Head Fixation Screws**
 - C. **Installation of the Upper Head Fixation Screws**
 - D. **SMS Flex Coil Installation (Optional)**
 - E. **Positioning the Patient**
 - F. **Removal of Patient**
 - G. **Removal from MRI Unit**
 - H. **Head Fixation Frame Cleaning**
 - I. **Cold Sterilization; Head Fixation Screws and Thumbwheels**

1. **Intended Use**

The Multi-Positional Head Fixation Frame is intended for use as a device to clamp and hold the patient's head in a particular position for procedures requiring Magnetic Resonance Imaging (MRI) of the brain structure and targets.

Warning: This device is intended for multiple uses. Contents are non-sterile.

Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

Note: The Head Fixation Frame is **MRI Conditional**

2. **Device Description**

These instructions for Use cover all ClearPoint Neuro Multi-Positional Head Fixation Frames.

The Multi-Positional Head Fixation Frame is designed to immobilize the patient's head during the surgical procedure. The patient's head is rigidly attached to the Head Fixation Frame by the physician in the desired position (lateral, oblique, prone or supine) for the planned surgical procedure.

See the accompanying scanner table specific Installation IFU for the frame provided which includes the contents of the assembly, associated devices, and instructions for attachment to the scanner table.

3. **General Warnings and Precautions**

The Head Fixation Frame is intended for multi-use and is provided non-sterile. It is to be used in a non-sterile and disinfected condition.

Warning: Do not clamp/affix the Head Fixation Frame to damaged or diseased bone. Only attach to stable bone to insure proper fixation.

Warning: Do not clamp/affix the Head Fixation Frame to patients with skull prostheses or previous reconstructive surgeries.

Warning: Do not over tighten the head fixation screws; select the proper clamping pressure related to the skull thickness and bone structure.

Warning: Make sure that the patient's eyes do not contact the skull pins during head positioning. Eye contact with the pins can lead to serious injury.

Caution: This device is designed to use the DORO Skull Pins, failing to use the recommend Skull Pins may result in improper function of the Head Fixation Frame. (Recommended: DORO[®] Skull Pins, Manufactured by PMI - Pro Med Instruments)

Caution: This device is to be used only by physicians trained by ClearPoint Neuro personnel.

Caution: The compatibility of other coil devices should be evaluated before use with the ClearPoint[®] System.

- Caution:** The Head Fixation Frame is designed to be robust; however, any attached Flex Coil (optional) is sensitive electronic equipment that requires handling with care.
- Caution:** All other coils must be removed from the MRI table, prior to installing the ClearPoint Neuro Head Fixation Frame.
- Caution:** Notify ClearPoint Neuro immediately if any screws break or holes are stripped out on the Head Fixation Frame.
- Note:** Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and appropriate Competent Authority of the Member State in which the incident occurred.
- Note:** Safe disposal of the device: The Titanium Skull Pins shall be treated as biohazardous materials and shall be disposed of accordingly. The Head Fixation Frame shall be disposed of per hospital policy.

General Precautions

Clean the Head Fixation Frame using standard hospital practices regarding cleaning, disinfecting and cold sterilization.

4. Use Instructions

A. Preparation

The Head Fixation Frame is provided non-sterile. The Head Fixation Frame uses Head Fixation Screws and Skull pins for patient fixation. The Skull Pins required to fixate the patient are single-use, sterile components, manufactured by:

Pro Med Instruments

(DORO® Skull Pins Titanium Tip: Item # 3006-20)

(DORO® Skull Pins Titanium Tip Children: Item # 3006-30)

WARNING: Do not use the Head Fixation Frame if any of its components are damaged.

1. Inspect the Head Fixation Frame along with the (optional) Flex Coil for defects. Do not proceed if defected items cannot be replaced.
2. Remove all accessory MRI coils from the MRI table (i.e. Spine Coil, Head Coil, and others).

B. Installation of the Lower Head Fixation Screws

Note: This process is common to all ClearPoint Neuro Multi-Positional Head Fixation Frames.

1. The Multi-positional Head Frame comes with two sets of Lower Fixation screws. There are short screws and long screws. The long screws are installed with the Head Frame. The short screws can be found in the clear plastic box that is included with the Head Frame.
2. Instructions for placing Lower screws in the Frame are given below.
3. Place Thumbwheel in the slot of the Head Fixation Frame. **(Figure: 1)**
4. Insert the squared end of the Lower Fixation Screw in the hole in the lower portion of the arc of the Head Fixation Frame. **(Figure: 2)**
5. Align the squared portion of the Thumbwheel with the squared end of the Lower Fixation Screw. **(Figure: 1)**


Figure 1

6. Retract the Lower Fixation Screw by rotating the thumb wheel. Advance the screw by rotating the thumb wheel in the opposite direction. **(Figure: 2)**


Figure 2

C. Installation of the Upper Head Fixation Screws

Note: This process is common to all ClearPoint Neuro Head Fixation Frames:

Insert the Upper Fixation Screws into the desired positions in the Head Fixation Frame. There are six available positions for the upper screws. Advance the screws until the ends are approximately 1" from the inside of the Head Fixation Frame (**Figure 3**).


Figure: 3

D. SMS Flex Coil Installation (Optional)

Place the SMS Flex Coil onto the upper and lower Head Fixation screws using the openings in the Flex Coil. (**Figure: 4**)


Figure: 4

Mount the four (4) disposable Skull Pins by removing them from the sterile package and placing them in the recessed ends of the Upper and Lower Fixation Screws. (Figure: 5)


Figure: 5

Recommend Skull Pins:

NGS-SP-01: Skull Pin Tips

MR Compatible Skull Pin Tips (DORO Skull Pins Titanium Tip: Item # 3006-20, Manufacturer: Pro Med)

NGS-SP-02: Pediatric Skull Pin Tips

MR Compatible Skull Pin Tips, Pediatric (DORO Skull Pins Titanium Tip Children: Item # 3006-30, Manufacturer: Pro Med)

Note: Handle the Skull Pins by the protection tube while removing from the package and placing into the Fixation Screws.

Note: Ensure each Skull Pin is fully seated.

E. Positioning the Patient

1. After fully seating the Skull Pins, remove the protection tubing from the Skull Pins and ensure the Fixation Screws and Skull Pins are retracted enough to be able to place the patients head into the Multi-Positional Head Fixation Frame.
2. Lay the patient on the MRI table with their head located in the Head Fixation Frame in the appropriate position, as determined by the Physician.
3. Raise the Lower Fixation Screws by rotating the Thumbwheels, counterclockwise to lift the patients head to the desired height in alignment with the Upper Fixation Screws and Skull Pins. **(Figure: 6)**
Warning: Make sure that the patient's eyes do not contact the skull pins during head positioning. Eye contact with the pins can lead to serious injury.

NOTE: For additional patient safety, eye protection during positioning may be utilized. Safety eye wear is available in both adult and pediatric sizes.


Figure: 6

4. Push firmly on the patient's Head to embed the bottom Skull Pins into the tissue.
5. Tighten the Upper Fixation Screws by rotating clockwise, until the appropriate force is applied to the patient's head. **(Figure: 7)**


Figure: 7

6. Tighten the Lower Fixation Screws by rotating the Thumbwheels counterclockwise until the appropriate force is applied to the patient's head.

Warning: Do Not Over tighten the Upper and Lower Fixation Screws.

Caution: The patient's head must remain immobile throughout the procedure.

F. Removal of Patient

Loosen the Upper Head Fixation Screws

1. Retract the Lower Fixation Screws.
2. Remove Patient from Head Fixation Frame.

G. Removal from MRI Unit

1. Remove the Skull Pins from the Head Fixation Frame and properly dispose.
Loosen the clamps that secure the head fixation frame to the table and carefully remove.

H. Head Fixation Frame Cleaning

The Head Fixation Frame is provided non-sterile and is designed for routine disinfecting. The Upper and Lower Fixation Screws are designed for multiple uses and routine disinfecting.

1. Remove the Upper and Lower Fixation Screws, along with the Lower Thumb Wheels. **(Figure: 8)**


Figure: 8

2. The entire system should be wiped down with appropriate cleaners.
3. Hardware should be soaked and cleaned with appropriate cleaners.
4. Flush all channels and threaded inserts.
5. Use a swab or similar to clean the four (4) threaded inserts within the Head Fixation Frame.
6. Dry the Head Fixation Frame with a lint free towel or equivalent.

I. Cold Sterilization; Head Fixation Screws and Thumbwheels

1. Clean the Head Fixation Screws and Thumbwheel as noted above.
 2. Activate the Cold Sterilization Solution per the manufacturer's instructions and Institutional policies.
 3. Immerse the Head Fixation Screws and Thumbwheels completely in the Cold Sterilization Solution.
 4. Agitate each Head Fixation Screw and Thumbwheel as it is immersed in the solution to eliminate air bubbles on the surface of the Head Fixation Screw.
 5. Ensure the solution fills the openings in the Head Fixation Screws.
 6. Leave the Head Fixation Screws and Thumbwheels immersed for the time recommended by the Cold Sterilization Solution manufacturer to obtain sterility.
 7. Using standard hospital sterile techniques remove the components from the solution and rinse thoroughly by immersing the components completely in three separate copious volumes of sterile saline or water. Each rinse should last a minimum of one minute. Use fresh sterile saline for each rinse.
- Caution:** Please refer to the Cold Sterilization Manufacturer's instructions, for complete instructions concerning the use of the solution as a sterilizing agent.

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ClearPoint Neuro, Inc.