CLEARPOINT

NEURO

SMARTGrid MR Planning Grid - BR

INSTRUCTIONS FOR USE

Table of Contents

I.	Intended Use
II.	Device Description
III.	General Warnings and Precautions
IV.	MRI Safety Information
V.	Use Instructions
A.	Preparation
B.	Marking the Burr Hole Site

- I. **Intended Use:** The SMARTGrid is an accessory to the ClearPoint System to assist in providing stereotactic planning, guidance and delivery of devices within the brain under Magnetic Resonance Imaging (MRI).

II. Device Description:

Package Contents:

NGS-SG-01-11-BR: SMARTGrid MR Planning Grid

Marking Grid (1)

Marking Tool (1)

This is a MR visible planning tool for surgical point of entry. The SMARTGrid Marking Tool is MR Conditional.

SMARTGrid

Marking Tool

III. General Warnings and Precautions

- Warning:** This device is intended for single patient use only. Contents of unopened, undamaged package are sterile. Do not re-sterilize.
- Warning:** There are no known and reliable means of cleaning, disinfecting, repairing, and sterilizing these devices that returns them to original specifications and renders them safe and effective for reuse.
- Caution:** Federal Law restricts this device to sale by or on the order of a physician.

IV. MRI Safety Information

MRI Safety Information: MR Conditional

A patient may be safely scanned under the following conditions. Failure to follow these conditions may result in an injury. Additional MRI Safety Information may be found at <https://www.clearpointneuro.com/eifus/> or by calling 888-287-9109.

Parameter	Condition of Use / Information
Static Magnetic Field Strength (B_0)	1.5T, 3T
Maximum Spatial Field Gradient (SPG)	22 T/m and 2,200 gauss/cm
RF Polarization	Circularly Polarized (CP)
RF Transmit Coil Type	Any Transmit RF Coil may be used
MR System (RF) Operating Modes or Constraints	Normal Mode
Whole-Body Averaged SAR	Whole-Body Averaged SAR ≤ 1.8 W/kg
Scan Duration and Wait Time	Scan for 15 minutes of continuous RF exposure with one or more MR imaging pulse sequences (scans or series) with 30 mins cooling period.
Artifacts	The presence of the implant may produce an image artifact.
Additional Information	The health state of the patient or the presence of other implants may require a reduction of the MRI limits. Applicable Components: SmartGrid Planning Grid (Marking Grid)

MRI Safety Information: MR Conditional

The device may be used safely in the MR environment under the following conditions. Failure to follow these conditions may result in an injury. Additional MRI Safety Information may be found at <https://www.clearpointneuro.com/eifus/> or by calling 888-287-9109.

Parameter	Condition of Use / Information
Static Magnetic Field Strength (B_0)	1.5T, 3T
Maximum Static Magnetic Field [mT] and [Gauss]	Do not exceed 158 mT (or 1580 Gauss) This corresponds to the distance of 639 mm from the MRI Scanner bore.
Additional Information	The health state of the patient or the presence of other implants may require a reduction of the MRI limits. Applicable Components: SmartGrid Planning Grid (Marking Tool)

Projectile Hazard
Equipment operation may be impacted
Do not exceed 158 mT (1580 Gauss)

V. Use Instructions:

A. Preparation

1. Position and secure patient in an appropriate head fixation frame and appropriate imaging coils to achieve desired image quality and provide access to the procedure site.
2. Install a sterile field MR Neuro Procedure Drape and attach the incise tape portion of the drape to the patient's head.


B. Marking Burr Hole Location

1. Remove paper back from Marking Grid.
2. Affix sticky side of Marking Grid onto the patient's head (over the incise tape portion of the drape), centered at the approximate location for the burr-hole.
3. Determine location on SmartGrid to mark skull per instructions from CLEARPOINT workstation (alpha-numeric coordinates).
4. Peel off top fluid filled layer of SmartGrid while firmly holding the lower layer in position.
5. Using Marking Tool and the coordinates, mark location on skull, entering through one of the holes provided on the Grid. The Marking Tool should be twisted/rotated several times using sufficient pressure to get a recognizable mark on the skull.

Caution: Only use minimum force required to leave a reference mark on the skull. This is only meant to provide surgeon a landmark, not to penetrate patient's skull.

Peel off remaining bottom layer of SMARTGrid.

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Beekley Corporation
ClearPoint Neuro, Inc.