SMARTGrid Planning Grid

INSTRUCTIONS FOR USE

EN

Table of Contents

- I. **Intended Use / Intended Purpose**
- II. **Device Description**
- III. **General Warnings and Precautions**
- IV. **MRI Safety Information**
- V. **Use Instructions**
 - A. **Preparation**
 - B. **Marking the Burr Hole Site**

- I. **Intended Use / Intended Purpose:** The SMARTGrid is an accessory to the ClearPoint System to assist in providing stereotactic planning, guidance and delivery of devices within the brain. The ClearPoint System is intended to provide stereotactic guidance for the placement and operation of instruments or devices during planning and operation of neurological procedures within an operating room environment and in conjunction with MR and/or CT imaging. During planning, the system is intended to provide functionality for the automatic identification, labeling, visualization, and quantification of segmentable brain structures from a set of loaded MR images. The ClearPoint System is intended as an integral part of procedures that have traditionally used stereotactic methodology. These procedures include biopsies, catheter and electrode insertion including deep brain stimulation (DBS) (asleep or awake) lead placement. When used in an MRI environment, the system is intended for use only with 1.5 and 3.0 Tesla MRI scanners and MR Conditional implants and devices.

Contraindications: The ClearPoint System is contraindicated for use with higher than 3.0 Tesla MRI scanners.

Contraindications for the use of the ClearPoint System and Accessory Kits are those contraindications generally accepted for intracranial surgery, including, but not limited to, severe, uncorrected bleeding disorder, local or systemic infection, and/or patient inability to tolerate general anesthesia.

Intended Users: The intended users are neurosurgeons who are trained in the use of these devices by ClearPoint Neuro personnel.

Intended Patient Population: Target populations for the ClearPoint System and Accessory Kits include any patient who is a candidate for precise and accurate stereotactic, image-guided placement of an intracranial device for therapeutic or diagnostic purposes.

II. Device Description:

Package Contents:

NGS-SG-01-11: SMARTGrid Planning Grid

Marking Grid (1), Marking Tool (1)

This is a MR visible planning tool for surgical point of entry. The SMARTGrid Marking Tool is MR Conditional.

SMARTGrid

Marking Tool
Herramienta de marcado
Outil de marquage
Strumento di
Markierungs-Tool
Narzędzie do znakowania
Mærkningsværktøj
Markeergereedschap
Markeringsverktyg
Ferramenta de Marcação

III. General Warnings and Precautions

- a. **Warning:** This device is intended for single patient use only. Contents of unopened, undamaged package are sterile. Do not re-sterilize.
- b. **Warning:** There are no known and reliable means of cleaning, disinfecting, repairing, and sterilizing these devices that returns them to original specifications and renders them safe and effective for reuse.

- c. **Warning:** Do not use if the package or any of the parts of the device are damaged.
- d. **Caution:** Federal Law restricts this device to sale by or on the order of a physician.

Note: Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and appropriate Competent Authority of the Member State in which the incident occurred.

IV. MRI Safety Information

MRI Safety Information: MR Conditional

A patient may be safely scanned under the following conditions. Failure to follow these conditions may result in an injury. Additional MRI Safety Information may be found at <https://www.clearpointneuro.com/eifus/> or by calling 888-287-9109.

Parameter	Condition of Use / Information
Static Magnetic Field Strength (B_0)	1.5T, 3T
Maximum Spatial Field Gradient (SPG)	22 T/m and 2,200 gauss/cm
RF Polarization	Circularly Polarized (CP)
RF Transmit Coil Type	Any Transmit RF Coil may be used
MR System (RF) Operating Modes or Constraints	Normal Mode
Whole-Body Averaged SAR	Whole-Body Averaged SAR ≤ 1.8 W/kg
Scan Duration and Wait Time	Scan for 15 minutes of continuous RF exposure with one or more MR imaging pulse sequences (scans or series) with 30 mins cooling period.
Artifacts	The presence of the implant may produce an image artifact.
Additional Information	The health state of the patient or the presence of other implants may require a reduction of the

	MRI limits. Applicable Components: SmartGrid Planning Grid (Marking Grid)
--	---

MRI Safety Information: MR Conditional

The device may be used safely in the MR environment under the following conditions. Failure to follow these conditions may result in an injury. Additional MRI Safety Information may be found at <https://www.clearpointneuro.com/eifus/> or by calling 888-287-9109.

Parameter	Condition of Use / Information
Static Magnetic Field Strength (B_0)	1.5T, 3T
Maximum Static Magnetic Field [mT] and [Gauss]	Do not exceed 158 mT (or 1580 Gauss) This corresponds to the distance of 639 mm from the MRI Scanner bore.
Additional Information	The health state of the patient or the presence of other implants may require a reduction of the MRI limits. Applicable Components: SmartGrid Planning Grid (Marking Tool)

Projectile Hazard
Equipment operation may be impacted
Do not exceed 158 mT (1580 Gauss)

V. Use Instructions:

A. Preparation

1. Position and secure patient in an appropriate head fixation frame and appropriate imaging coils to achieve desired image quality and provide access to the procedure site.
2. Install a sterile field drape and attach the incise tape portion of the drape to the patient's head.

B. Marking Burr Hole Location

1. Remove paper back from Marking Grid.
2. Affix sticky side of Marking Grid onto the patient's head (over the incise tape portion of the drape), centered at the approximate location for the burr-hole.
3. Determine location on SmartGrid to mark skull per instructions from CLEARPOINT workstation (alpha-numeric coordinates).
4. Peel off top fluid-filled layer of SmartGrid while firmly holding the lower layer in position.
5. Using Marking Tool and the coordinates, mark location on skull, entering through one of the holes provided on the Grid. The Marking Tool should be twisted/rotated several times using sufficient pressure to get a recognizable mark on the skull.

Caution: Ensure the top fluid-filled layer is removed prior to using the Marking Tool. Do not use the Marking Tool to mark the location through the fluid-filled layer.

Caution: Only use minimum force required to leave a reference mark on the skull. This is only meant to provide surgeon a landmark, not to penetrate patient's skull.

6. Peel off remaining bottom layer of SMARTGrid.

Caution: Ensure all pieces of the grid layer are fully removed from the patient's head.

Note: Safe disposal of the device: The device shall be treated as biohazardous materials and shall be disposed of accordingly per hospital policy.

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Beekley Corporation
ClearPoint Neuro, Inc.